

Salivary Melatonin ELISA Kit

Item #: 1-3402 (single), 1-3402-5 (5-pack)




The Most **Accurate, Easiest** and **Fastest** Melatonin Assay!

Assay Comparison:		Salimetrics Assay	Other Assays
Accurate:	Antibody Designed for Results in Saliva	✓	-
Fast:	Results in < 4 Hours	✓	-
Accurate:	Low Analytical Sensitivity	✓	✓
Fast:	No Sample Pre-Treatment (One Wash Step)	✓	-
Accurate:	Consistent, Reliable Results	✓	-
Easy:	100 µL Sample Volume	✓	-
Accurate:	Average %CV < 5%	✓	-
Easy:	Simple Protocol	✓	-
Accurate:	Best Precision for DLMO Levels	✓	-

Protocol Comparison

Assay Procedure	Time Required (hrs)	
	Other	Salimetrics
Sample Extraction, Wash, Diluent	0.4	-
Pipette Samples, Stds, Controls	0.2	0.2
Incubation #1 & Wash	20.0	3.0
Conjugate, Incubation #2 & Wash	3.1	-
SA-HRP, Incubation #3 & Wash	1.1	-
TMB, Incubate, Stop, Read	0.6	0.6
Total Time (hrs)	25	4

Trust Your Assay Always Accurate DLMO Determination



* Samples run in triplicate



Better Collection. Better Assays. Better Results. Better Science.